

INVACARE® MATRIX®

# PSP (Posture Seat Polymer) CUSHION

---



The **Invacare Matrix PSP** (Posture Seat Polymer) cushion combines the stability and support of the PS foam base with a polymer overlay for extra protection and shear reduction. An ultra-lightweight cushion providing zero maintenance, worry free skin protection and positioning.

## Call to order

**1.800.668.5324**  
invacare.ca

**1.866.748.7943**  
motionconcepts.com



**matrix**  
SEATING SERIES



**Features and Functional Benefits**

- Startex™ outer cover – Reversible can be used smooth or fabric side out
- Moisture resistant and breathable to keep cool and dry
- Includes hook and loop strips
- Startex™ inner cover for extra moisture protection
- Molded high resilient polyurethane foam made with Ultra-fresh™ anti-microbial
- Rear cushion radius and waterfall front edge
- Gently contoured to provide neutral leg position
- Molded polymer insert for extra protection and shear reduction
- Unique waffled relief zone allows immersion in ischial area
- Designed to provide off-loading to hips and thighs for support and stability

ADP# SE0001583

**Optional Accessories**



**POK – Pelvic Obliquity Kit**  
(2 pieces) accommodates up to one inch of obliquity



**Cushion Rigidizer –** lightweight and low profile

**Invacare Corporation**  
www.invacare.com

**USA**  
One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900


**Canada**  
570 Matheson Blvd. E.  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

**Motion Concepts**  
www.motionconcepts.com  
info@motionconcepts.com

**USA**  
700 Ensminger Rd., Suite 112  
Tonawanda, NY 14150  
(888) 433-6818  
Fax: (888) 433-6834

**Canada**  
84 Citation Dr., Unit 1-7  
Concord, ON L4K 3C1  
(866) 748-7943  
(905)-695-0134  
Fax: 905-695-0138

© 2018  
Invacare Corporation.  
All rights reserved.  
Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Matrix is a registered trademark of Motion Concepts. Otto Bock is a trademark of Otto Bock HealthCare. Specifications are subject to change without notification.  
Form No. 15-034  
Rev. 03/19



Model #	Width	Depth	Height Abductor/Adductor	Height at Leg Trough	Product Weight	Weight Capacity	Limited Warranty
PSP (Standard)	14", 16", 17", 18", 20"	16", 18", 20", 22"	3.5"/3.5"	2.75"	2.2 lb. and up	300 lb.	3 Years

