



Invacare® Scooters
Suggested Canadian Price List
Effective March 2nd, 2009
Updated February 1st, 2012
***UPDATED - DENOTES CHANGES**

NAME _____ PHONE NO _____

ACCOUNT # _____ P.O. _____ DATE _____

ADDRESS _____

INVACARE CANADA SCOOTERS

LEO™ SCOOTERS

- LEO-4S 4-Wheel Midsize, Silver \$2,999
- LEO-4B 4-Wheel Midsize, Blue \$2,999

LYNX™ SCOOTERS

- L-3R 3-Wheel Microportable, Red \$1,441
- ~~L-3B~~ ~~3-Wheel Microportable, Blue - UPDATED - 1/15/2012 - NLA~~ ~~\$1,441~~
- L-4R 4-Wheel Microportable, Red \$1,544
- L-4B 4-Wheel Microportable, Blue \$1,544
- ~~L-3XR~~ ~~3-Wheel Compact, Red - UPDATED - 1/15/2012 - NLA~~ ~~\$1,956~~
- ~~L-3XB~~ ~~3-Wheel Compact, Blue - UPDATED - 1/15/2012 - NLA~~ ~~\$1,956~~

PEGASUS

- PEGASUS3S Pegasus 3 Wheel Scooter with 5A Charger, Silver \$3,699
- PEGASUS3B Pegasus 3 Wheel Scooter with 5A Charger, Blue \$3,699
- PEGASUS4S Pegasus 4 Wheel Scooter with 5A Charger, Silver \$3,799
- PEGASUS4B Pegasus 4 Wheel Scooter with 5A Charger, Blue \$3,799

COMET

- COMET4S Comet 4 Wheel Scooter with 8A Charger, Silver \$4,299
- COMET4B Comet 4 Wheel Scooter with 8A Charger, Blue \$4,299

COMET

- COMET-HD Comet 4 Wheel Heavy Duty Scooter w/ 8A Charger, Grey \$4,899

Leo Features and Options
Choice of Silver or Blue Finishes Flat Free Tires Built-In Anti-Splash Guards Infinitely Adjustable Tiller 36AH Batteries Included
Lynx Features and Options
Choice of Red or Blue Finishes Horn Flat Free Tires Self Diagnostics On/Off Board Charger on L-3 / L-4 Batteries Included
Pegasus Features and Options
Choice of Silver or Blue Finishes Adjustable rear suspension. Front spring suspension 4-Pole Motor Automatic Slow Down Technology Built-In Anti-Splash Guards Infinitely Adjustable Tiller Pneumatic Tires
Comet/Comet HD Features and Options
Choice of Silver or Blue (Comet), Slate Grey (HD) 4-Wheel independent Adjustable Suspension Automatic Slow Down Technology Built-In Anti-Splash Guards Infinitely Adjustable Tiller 15 km/h (9.5 mph speed) Comet 10 km/h (6 mph speed) Comet HD Pneumatic Tires

	LYNX		LEO	PEGASUS		COMET	
	3-Wheel Micro-Portable	4-Wheel Micro-Portable	4-Wheel Midsize	3-Wheel	4-Wheel	4-Wheel	4-Wheel Heavy Duty
Overall Dimensions							
Width	20"	20"	22"	24 1/2"	26"	26"	26" - 30" (inc. seat)
Length	39.8	39.8	47"	50"	51"	58"	58"
Height	36"	36"	36"	47 1/2"	47 1/2"	49 1/2"	58"
Step Height	4"	4"	3"	4"	4"	5"	5"
Usable Foot Space	15 1/4"	15 1/4"	12"	11 1/2"	12 1/2"	13"	13"
Seat Dimensions							
Seat-to-Floor	14" - 16"	14" - 16"	18" - 22"	16 1/2" - 19 1/2"	16 1/2" - 19 1/2"	16 1/2" - 19 1/2"	19" - 21"
Width	16 1/2"	16 1/2"	19"	20"	20"	20"	26"
Depth	14 1/2"	14 1/2"	16"	18 1/2"	18 1/2"	18 1/2"	22 1/2"
Back Height w/o Headrest	16"	16"	18 1/2"	20"	20"	20"	32"
Back Height w/Headrest	N/A	N/A	N/A	27" - 29 1/2"	27" - 29 1/2"	27" - 29 1/2"	32" - 35"
Back Angle	N/A	N/A	N/A	95° - 145°	95° - 145°	95° - 145°	95° - 120°
Arm Dimensions							
Length	11.6"	11.6"	14"	14"	14"	14"	14"
Height Between Seat and Arms	9"	9"	7"	8" - 10"	8" - 10"	8" - 10"	5" - 7"
Width Between Arms	15" - 18"	15" - 18"	23" - 27"	20" - 23"	20" - 23"	20" - 23"	26" - 29"
Battery	2 x 12 Ah (inc.) On or Off Board	2 x 12 Ah (inc.) On or Off Board	2 x 36 Ah (inc.) Off Board, 24V, 5 Ah	40 Ah, 2 x 12 V Off Board, 24V, 5 Ah	40 Ah, 2 x 12 V Off Board, 24V, 5 Ah	70 Ah, 2 x 12 V Off Board, 24 V 8 Ah	70 Ah, 2 x 12 V Off Board, 24 V 8 Ah
Charging System	On or Off Board	On or Off Board	Off Board, 24V, 5 Ah	Off Board, 24V, 5 Ah	Off Board, 24V, 5 Ah	Off Board, 24 V 8 Ah	Off Board, 24 V 8 Ah
Speed	8 km/h (5 mph)	8 km/h (5 mph)	8 km/h (5 mph)	12 km/h (7.5 miles)	12 km/h (7.5 mph)	15 km/h (9.5 mph)	10 km/h (6 mph)
Estimated Travelling Range*	11 km (7 miles)	11 km (7 miles)	36 km (22 miles)	42 km (26 miles)	42 km (26 miles)	55 km (34 miles)	50 km (31 miles)
Weight Capacity (lbs.)	300	300	350**	330	330	390	485
Maximum Incline Capability	8°	8°	8°	10°	12°	10°	10°
Suspension	N/A	N/A	N/A	Rear independent adjustable, front spring Standard	Rear independent adjustable, front spring Standard	Four wheel independent adjustable Standard	HD four wheel independent adjustable Standard
Traction Control System	N/A	N/A	N/A	Magnetic and dual front manual	Magnetic and dual front manual	Magnetic and dual front manual	Magnetic and dual front manual
Braking	Disc	Disc	Disc	Programmable	Programmable	Programmable	Programmable
Electronics	Programmable	Programmable	Programmable	Programmable	Programmable	Programmable	Programmable
Motor	2-Pole, 150 W	2-Pole, 150 W	4-Pole, 240 W	4-Pole, 250 W	4-Pole, 250 W	4-Pole, 250 W	4-Pole, 650 W
Tires	Foam filled	Foam filled	Foam filled	Pneumatic	Pneumatic	Pneumatic	Pneumatic
Front	8" x 2"	8" x 2"	9" x 2"	11" x 3 1/2"	11" x 3 1/2"	13" x 4"	13" x 4"
Rear	8" x 2"	8" x 2"	9" x 2"	11" x 3 1/2"	11" x 3 1/2"	13" x 4"	13" x 4"
Turning Radius	39.4"	51.2"	51.2"	44"	55"	64"	64"
Weight w/o batteries (lbs.)	103	112	184	156	172	185	212
Warranty	Two year limited warranty against manufacturing defects (excluding upholstered materials, padded materials and tires/wheels). Six months on batteries (Lynx, Leo). Manufacturer's warranty does not cover normal wear.						

