



NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_

ACCOUNT # \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_

ADDRESS \_\_\_\_\_

**STOCK WHEELCHAIR CONFIGURATIONS**

- TREX26R 16" X 16" Frame with Removable Fixed Height, Desk Length Arms \$612
- TREX26RF 16" X 16" Frame with Removable Fixed Height, Full Length Arms \$612
- TREX28FF 18" X 16" Frame with Permanent Arms & Fixed Footrests \$569
- TREX28P 18" X 16" Frame with Permanent Arms \$528
- TREX28R 18" X 16" Frame with Removable Fixed Height, Desk Length Arms \$612
- TREX28RF 18" X 16" Frame with Removable Fixed Height, Full Length Arms \$612
- TREX20P 20" X 16" Frame with Permanent Arms \$528
- TREX20R 20" X 16" Frame with Removable Fixed Height, Desk Length Arms \$612
- TREX20RF 20" X 16" Frame with Removable Fixed Height, Full Length Arms \$612

**STOCK FRONT RIGGINGS**

- T93HE Economy Hemi Footrests; Composite Footplates without Heel Loops \$111
- T93HC Hemi Footrests, Composite Footplates with heel loops \$152
- T93HA Hemi Footrests, Aluminum Footplates with heel loops \$152
- T94HE Economy Hemi Elevating Legrests; Composite Footplates & Calf Pads \$232
- T94HC Hemi Elevating Legrests; Composite Footplates & Padded Calf Pads \$290
- T94HA Hemi Elevating Legrests; Aluminum Footplates & Padded Calf Pads \$290
- AHL4A Hemi Smartleg Articulating Legrests; Aluminum Footplates & Padded Calf Pad \$377

**STOCK OPTIONS & ACCESSORIES**

**ARMS**

- 19ST Full Length Adjustable Height Arms (10" – 14" Adjustable Range) \$107
- 29ST Desk Length Adjustable Height Arms (10" – 14" Adjustable Range) \$107

**ANTI-TIPPERS**

- 9758ST Rear Anti-Tippers for TREX28FF \$117
- 1360 Adjustable Anti-Tippers (Does not fit fixed frame model) \$105

**MISCELLANEOUS**

- 1329BKST 8" Wheel Lock Extensions \$61
- 1335ST Non-Folding Device \$61
- 1514 Overhead Anti-Theft Device \$72
- 1515ST Push Button Style Seat Positioning Strap \$81
- 1496 IV/Oxygen Cylinder Holder Package Assembly \$254 **Footnotes:**
- 1497 Oxygen Cylinder Holder Package Assembly \$190 1 Not available on Fixed Frame model.
- 1820ST Non-Removable Device for Front Riggings \$39 2 Not available on 20" widths.

All Stock Tracer EX2 Wheelchairs Include: Standard	
* Dual Axle Positions (17 ½" and 19 ½" Seat-to-Floor)	
* 16" Back Height	
* 8" x 1" Solid Rubber Casters	
* 24" Composite Wheels with Urethane Tires	
* Composite Handrims	
* Midnight Blue Vinyl Upholstery (U550)	
* Padded Armpads	
* Wheel Locks	

Specifications	
Seat Widths	16", 18", or 20"
Overall Width	Permanent Arm – Seat Width plus 7 ¼" Removable Arm – Seat Width plus 8 ½"
Seat Depth	16"
Overall Depth	31 7/8" (18x16 without Front Riggings)
Back Height	16"
Seat-to-floor	17 ½" or 19 ½" (21" on Fixed Frame)
Weight Limit	250 lbs.
Chair Weight	36 lbs. (18x16 without Front Riggings)
Shipping Weight	45 ½ lbs.

**Invacare Canada reserves the right to change prices and specifications without notice.**

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324(Canada) Fax: (905) 890-5244 (800) 668-5478(Canada)  
 Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate.  
 For the most current pricing information visit www.invacare.ca.