

Top End® TRANSFORMER™ ALL Sport

Suggested Canadian Price List

Effective March 1st, 2005 Revised January 6th, 2009
 *UPDATED – DENOTES CHANGES



For accuracy of ordering, please fax order to Invacare Canada LP @ 905-590-5244 or
 Toll Free @ 800-668-5478

www.topendwheelchair.com

Provider Name _____ Date _____ Top End Job Number _____
 Account # _____ P.O. # _____ Phone # _____ Fax # _____
 Address _____ City _____ Prov. _____ Postal Code _____
 Customer Name _____ Phone No _____
 Address _____ City _____ Prov. _____ Postal Code _____
Ship order to: Provider Customer **Order will ship to Provider unless otherwise noted.**

CHAIR TYPE

- TRAF Transformer All Sport**3495.00
- Seat Width:** 12" 13" 14" 15" 16" 17" 18" 19" 20"
- Seat Depth:** 12" 13" 14" 15" 16" 17" 18" 19" 20"
- Adjustable Frame Height: choose tall or short**
- TALL:** Front Seat Height 18-21" Seat at _____
 Rear Seat Height 15-21" Seat at _____¹
- SHORT:** Front Seat Height 16-18" Seat at _____
 Rear Seat Height 15-18" Seat at _____¹
- Back Type/Angle: (choose one)**
- Folding back / Adjustable angle
- Fixed back / 90° Angle
- Back Height:** 8-11" 10-14" 12-16" 14-18"
- Wing/bumper: (choose one)**
- Full wing (available with camber of 12°, 15°, 18°, or 20°)
- Bumper
- TRA190** Both wing and bumper165.00
- Back Upholstery** Fold Over**Standard**
- TRA49** Adjustable tension w/straps165.00
- Camber System**
- Single Camber (0-20) with camber tube**Standard**
 0° 3° 6° 9° 12° 15° 18° 20°
- TRA64** Extra camber tube230.00
 0° 3° 6° 9° 12° 15° 18° 20°
- TRA66** Adjustable camber system (choose 2)390.00
 Primary angle
 0° 3° 6° 9° 12°
 Secondary angle
 0° 3° 6° 9° 12°
- TRA46** Additional camber inserts (choose 2)165.00
 0° 3° 6° 9° 12°
- TRA48** Quick release lever40.00
- High performance wheels/options:Standard**
 24", 25", 26" and 700 C wheels are double wide rims,
 high flange hub equipped with high pressure clincher
 tires and quick release axles.
- Rear wheel size** 22" (25-501) 24" (25-540)
 25" (25-559) 26" (25-590) 700C
- TRA67** Spinery® Spox wheels upgrade695.00
 22" (25-520) 24" (25-540) 25" (25-559)
 26" (25-590) 700C
 Spoke color Black Yellow²
- TRA68** Kenda® tires Red Yellow90.00³
- TRA6** Crossed spokes upgrade165.00⁴
- TRA21** Spoke guards95.00²
 Black Clear
- TRA27** Quad release axles165.00
- TRA77** Titanium axles235.00

Pushrim type:

- Aluminum**Standard**
- TRA8** Coated Pushrims125.00
- TRA75** Titanium pushrim610.00⁵

Pushrim mount:

- Close mount **Standard Mount**

Front wheels:

- 3" Solid urethane w/ matching mini-fork**Standard**
- TRA15** Lighted front wheels and anti-tip165.00

Footrest:

- Adjustable height/angle one piece footplate**Standard**
- TRA111** Adj ht/angle one-piece footplate w/toe clips295.00⁶
- TRA191** Adaptive footplate spacers 2" 4" 6"130.00⁶
- Calf strap No Charge

Rear Anti-Tip

- Single swivel**Standard**
- TRA193** Double swivel165.00

SAFETY AND COMFORT OPTIONS:

- TRA14** Cushion: soft 2" 3" 4"100.00
- TRA33** Cushion: hard 2" 3" 4"125.00
- TRA194** Solid seat65.00
- TRA18** Removable sideguards165.00
- TRA51** Folding sideguards: Rigid with slot175.00⁷
- TRA121** Fold-in style sideguards170.00⁷
- TER95** Clamp-On Fixed aluminum side guards125.00
- TRA35** Frame protector (rear)30.00
- TRA37** Frame protector pair (front)55.00
- TRA94** Alignment gauge120.00
- TER126** Click strap large (pad length-12")105.00
- TER125** Click strap med (pad length-10")105.00
- TER124** Click strap small (pad length-8")105.00
- TRA17** Positioning straps – **Hook and Loop**60.00
- TRA19** Positioning straps – **Auto Style Buckle**75.00⁷
- TRA2** Wheel locks – **Push to Lock Hi-Mount**300.00⁷
- TRA3** Wheel locks – **Scissor Style Low-Mount**250.00⁸
- Race22** Top End T-Shirt L XL15.00 NET⁹
- TRA120** Fencing arm165.00

FRAME COLOR AND UPHOLSTERY

- UPHOLSTERY** Black std UPHI Camouflage230.00

CONTEMPORARY COLORS

- Wet Black Black Prism White Snow
- Silver Metallic Silver Vein Copper Penny

WAY COOL COLORS

- Electric Blue Midnight Blue Cosmic Blue
- Electric Teal Emerald Green Grasshopper Green

TOO HOT COLORS

FOOTNOTES

- ⁵ Titanium Pushrims not available w/700C or 22" wheels.
- ⁶ Needed for lower leg lengths of less than 15" and 12" (tall and short frames respectively).
- ⁷ Folding back only. Folding sideguard rigid w/slot not avail. With hi-mount wheellocks.
- ⁸ Scissor style low-mount not available with camber over 6 degrees. Other restrictions may apply depending on configuration.
- ⁹ Dealer discount does not apply.

FOOTNOTES

- ¹ Rear seat height is dependent on rear wheel size. 15" RSH requires 22" or 24" rear wheels.
- ² If no color is chosen, black will be shipped.
- ³ Kenda tires are available in 24", 25", and 26" sizes.
- ⁴ Crossed Spokes not available with Spinery wheels.

Top End® TRANSFORMER™ All Sport

- Tangerine Red Electric Red Grape Madness
 Safety Yellow Sunkissed Orange Fruit Fusion

TRANSFORMER OPTIONS

- | | |
|--------------------------------------|------------------------|
| * Removable Sideguards | * Spokeguard |
| * Folding Sideguards | * Spinergy Wheels |
| * Double Swivel Anti-Tippers | * Titanium Axles |
| * Clamp-on Raised Footrest | * Titanium Handrims |
| * Footplate Spacers | * Titanium Fork Stems |
| * Cushions | * Coated Handrims |
| * Solid Seat | * Frame Protectors |
| * Adjustable Tension Back Upholstery | * Extra Camber Inserts |

Special Instructions:

TRANSFORMER PRODUCT PARAMETERS

Adjustable 4 Wheel Multi-Use Sports Chair

Weight Limit	250 lbs
Material:	Aluminum Frame/Chrome-moly wing or bumper
Seat Widths:	12" – 20"
Seat Depths:	12" – 20"
Seat Upholstery:	Adjustable Tension
Back Upholstery:	Fold Over
Rear Wheels:	22-700c High Performance
Tires:	High Performance Clincher
Front Wheels:	3" Solid
Adjustable Rear Seat Height:	Tall or Short
Adjustable Front Seat Height:	Tall or Short
Back:	Folding Back Standard
Fixed Back:	No Charge Option
Back Height:	Adjustable
Back Angle:	Fixed Back is 90° to Floor Folding Back is Adjustable in 5° Increments
Center of Gravity:	Adjustable
Camber:	Camber Bar 2000 0° - 20°
Single Swivel Anti-Tip:	Adjustable
Height/Angle Footplate:	Standard
Adjustable Height/Angle Footplate:	Standard
Removable Full Wing or Bumper	Standard
Calf Strap	Standard
Paint:	All Basic Top End Colors
Frame Length is determined by Seat Depth plus 9" for Folding Back	
Frame Length is determined by Seat Depth plus 7" for Fixed Back	
Footrest Width is determined by Seat Width Minus 5"	

NOTE: All specifications and dimensions are approximate.

Invacare Canada LP

HEAD OFFICE –

570 Matheson Blvd, Unit #8,
Mississauga, ON
L4Z 4G4
Phone: (905) 890-8300 (800)
668-5324
Fax: (905) 890-5244 (800)
668-5478

Invacare Canada LP

Regional Centre - 16769

Boul. Hymus, Kirkland, P.Q.
H9H 3L4
Phones: (514) 630-6080 (800)
363-0137
Fax: (514) 630-6120

www.invacare.ca

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Spinergy is a registered trademark for Spinergy Corp.
Kenda is a registered trademark of Kenda Rubber Industrial Co. Ltd.
Specifications and prices are subject to change without notification.
©2007 Invacare Corporation